


## OVERVIEW OF THE CALIBRATION PROCESS

- A. Add and tighten all the objectives.
- B. Verify the camera to stage alignment.
- C. Calibrate lenses with calibration grids.**
- D. Perform parcentric/parfocal calibration.
- E. Check the accuracy of the calibration.




In this tutorial, we demonstrate how to calibrate 5x, 40x, and 63x lenses:


## CALIBRATE LENSES WITH CALIBRATION GRIDS

1. Load the calibration grid slide onto the microscope.
2. Select the lowest power objective (5x in this tutorial) on the microscope.
3. Enable **Joy Track** (click the **Joy Track** icon) to locate the 250  $\mu\text{m}$  (large) grid.
4. Turn off **Joy Track** (click the **Joy Track** icon again).
5. Verify that the software lens is set to 5x to match the objective.

*If there is no software lens defined for the objective, see [Calibration: Creating a Lens](#).*




6. Select **Tools>Grid Tune Current Lens** to open the **Grid Tune Current Lens** window.
  - ▶ Enter the box size that corresponds to the grid you're using: 250  $\mu\text{m}$  for 5x




7. A white grid with dotted lines appears. Align the white grid with the slide grid (black lines).  
*If the slide grid appears to be lopsided compared to the white grid on the screen, see [Troubleshooting](#).*

a. Drag the anchor to the corner of a cell located in the top left of the screen.




b. Align the vertical lines.

- i. Hover the mouse over a white vertical line until the cursor changes to a double arrow.
- ii. Drag the white dotted line to align it with the closest vertical slide grid line.
- iii. Repeat i-ii until all vertical lines are aligned.



c. Align the horizontal lines.

- i. Hover the mouse over a white horizontal line until the cursor changes to a double arrow.
- ii. Drag to align with the closest horizontal slide grid line.
- iii. Repeat i-ii until all horizontal lines are aligned.



- When the white grid is aligned as accurately as possible with the slide grid, right-click and select **Finish Calibrating Current Lens**.

*The 5x lens is calibrated.*



- Enable **Joy Track** to locate the 25  $\mu\text{m}$  (small) grid, then turn off **Joy Track**.
- Select the next objective (40x in this tutorial) on the microscope and verify that the software lens is set to 40x to match the objective.
- Enable **Joy Track** to move inside the 25  $\mu\text{m}$  grid toward the top left corner, then turn off **Joy Track**.
- Select **Tools>Grid Tune Current Lens** to open the **Grid Tune Current Lens** window.
  - Enter the box size that corresponds to the grid you're using: 25  $\mu\text{m}$  for 40x
- Another white grid with dotted lines appears. Align the white grid with the slide grid (black lines).
  - Drag the anchor to the corner of a cell located in the top left of the screen.
  - Align the vertical lines.
  - Align the horizontal lines.

- When the white grid is aligned as accurately as possible with the slide grid, right-click and select **Finish Calibrating Current Lens**.

*The 40x lens is calibrated.*



- Select the next objective (63x in this tutorial) on the microscope and verify that the software lens is set to 43x to match the objective.

*Don't use oil.*

- Select **Tools>Grid Tune Current Lens** to open the **Grid Tune Current Lens** window.
  - Enter the box size that corresponds to the grid you're using: 25  $\mu\text{m}$  for 63x

- Repeat steps 13-14.

Once the lenses are calibrated, you're ready to perform the parcentric/parfocal calibration.