

# NeuroInfo® software

Release Notes

## Version 2021.1.1 (November 9, 2020)

### NEW FEATURES AND ENHANCEMENTS

- Automated cell detection capabilities have been expanded
  - Improvements to the speed and accuracy
  - Now includes neural-network classifiers for intelligent identification and counting of specific neuronal subtypes
- A new, smart image background correction enables clear visualization and interpretation of compiled brain volumes
- This release employs the new high-performance imaging engine that reads, renders and writes big microscopy image and trace data at speeds up to 5000x faster than previous versions