



NeuroInfo®

Release Notes

VERSION 2023.2.1 (NOVEMBER 2023)

NEW FEATURES AND ENHANCEMENTS

- Added a new analysis tool for reporting experimental brain region volumes and intensities; these results can be exported to Excel
- In the Brain Region Ontology list, subregions can be automatically selected when the region is selected
- Markers can now be placed in the 3D visualization window
- Improved rendering and display responsiveness for millions of markers in 3D views
- Increased the speed of *Map Experiment to Atlas* for millions of markers from hours to minutes

ISSUES RESOLVED

- Improved support for image file formats
- Fixed the radius text field in the *Map Experiment to Atlas* tool
- Load and Save Transform options are now available after a new file is created