



Vesselucida® Explorer

Release Notes

Version 2021.1.1 (2/25/2021)

NEW FEATURES AND ENHANCEMENTS

- Puncta analyses added, including colocalization
- Easier software authorization: Vesselucida Explorer software autofills the online Authorization Request Form with information about your computer system. Just fill in a few fields with information about you and your lab to submit the form.
- Added ability to shell contours in the YZ and XZ orientation
- Units of measure were added to the Loops analysis
- Marker Detail analysis is now included for all markers
- Tube wall results in Contour analysis now include g-ratio computation
- New report on distance from individual markers to contours
- New report on distance from individual markers to vessel segments
- Added help links from Analyses window

ISSUES RESOLVED

- Fixed crash in Batch analysis
- Fixed issue that occurred when opening additional data files
- Fixed issue with the display of Closed Surface Analysis and contours
- Fixed issue with marquee select