



Vesselucida® Explorer

Release Notes

Version 2022.1.1 (December 2022)

NEW FEATURES AND ENHANCEMENTS

- New Sholl analysis for vessels
- Added help links from Analyses window
- Centroid coordinates added to the Tube wall analysis table
- Nearest Neighbor puncta analysis report expanded to all objects

ISSUES RESOLVED

- Fixed crash in Batch analysis
- Resolved issue using Tube Wall analysis with other analyses
- Puncta are highlighted when using Select all
- Resolved issue with calculating contours in closed contours
- Fixed issue with Closed Surface analysis working with the Batch export